

Silicone	Straight, Coiled	15F	Double
Material	Tip Designs	French Size	Cuff

COILED TIP
Provides a gentler, less traumatic, showering effect of dialysate on the peritoneum and prevents catheter migration.



DOUBLE CUFF
Ensures optimal outcomes by allowing one cuff to be placed in the rectus muscle and the other in the subcutaneous tissue.



PREFORMED BEND
The lateral or downward direction minimizes the risk of exit-site infection.



	CATHETER ONLY - STRAIGHT	
MC20VS43LC	STRAIGHT, 43CM LEFT V-NECK WITH DOUBLE CUFF	5/BOX
MC20VS43RC	STRAIGHT, 43CM RIGHT V-NECK WITH DOUBLE CUFF	5/BOX

	BASIC SET	
MC20VS43LS	STRAIGHT, 43CM LEFT V-NECK WITH DOUBLE CUFF	5/BOX
MC20VS43RS	STRAIGHT, 43CM RIGHT V-NECK WITH DOUBLE CUFF	5/BOX

	CATHETER ONLY - COILED	
MC20VC63LC	COILED, 63CM LEFT V-NECK WITH DOUBLE CUFF	5/BOX
MC20VC63RC	COILED, 63CM RIGHT V-NECK WITH DOUBLE CUFF	5/BOX

	BASIC SET	
MC20VC63LS	COILED, 63CM LEFT V-NECK WITH DOUBLE CUFF	5/BOX
MC20VC63RS	COILED, 63CM RIGHT V-NECK WITH DOUBLE CUFF	5/BOX

	ACCESSORIES	
MP217	17F VASCU-SHEATH INTRODUCER SET	5/BOX
MC50	MALE LUER LOCK END CAP	100/BOX
MC51	BETA CLAMP	100/BOX
MBCA	BARBED LUER LOCK CATHETER ADAPTOR	20/BOX
MCS	CATHETER STYLET	5/BOX
MTS	TUNNELING STYLET	10/BOX
40501	PREOPERATIVE MARKING STENCIL PACK	1/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter
- (1) Female Luer
- (1) Clamp
- (1) Injection Cap

BASIC SET CONTAINS:

- (1) Catheter
- (1) Female Luer
- (1) Clamp
- (1) Introducer Needle
- (1) Guidewire
- (1) Scalpel
- (1) Peelable Sheath Introducer
- (1) Tunneler
- (1) 10cc Syringe
- (1) Injection Cap
- (6) 4" x 4" Gauze

PD STENCIL PACK
IFU, Preoperative Mapping Record & Stencil

